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nal von minden Drug-Screen® Dip Test: Test Procedure & Evaluation

Test Procedure Single Test

- 1** To prepare you will need: test strips, urine cup with urine sample, gloves

- 2** Dip the test strips into the urine sample for 15-30 seconds, until the test zone is dyed a consistent pale yellow. Make sure that the urine level does not exceed the horizontal red line.

Test Procedure Multi-Dip Test

- 1** To prepare you will need: test strips, urine cup, urine sample, gloves.

- 2** Dip the test strips into the urine sample for 15-30 seconds, until the test zone is dyed a consistent pale yellow. Make sure that the urine level does not exceed the horizontal red line.

- 3** Lay the test strip on a clean and even surface

- 4** Read the results after 1-5 minutes. Do not evaluate results after 10 minutes, as they are no longer valid.

1-5 minutes

Evaluation of Results

Negative:

In the result window, two red lines will appear, one in the control area (C) and one in the test range (T).

The color strength of the test line (T) can vary. This result indicates that the concentration of the drug is under the detection limit, or that the drug is not present at all.

Positive:

In the result window, only one red line appears in the control area (C) and no line appears in the test range (T).

Invalid:

No control line (C) appears in the result window.